

PUSEN

mediform

## Single Use Cysto-Nephro Videoscope

PC200-AS | PC200-AR | PC200-S | PC200-R

*Always Sterilized  
Always Available  
Always Works*

Improved Clinical Outcomes



# PC200

## Straight Insertion Channel

Smooth navigation to low the risks of device damage



## Improved Surgeon Control

Manageable integrated suction button for optimal irrigation flow with superior intraoperative vision



High Quality Imaging



Suction Control



High Brightness



Straight Channel

Optical System	
Resolution	720 * 720 pixels
Field of View	90°
Depth of Field	5-80 mm
Insertion Portion	
Angle of Deflection	210° up/down
Diameter of Insertion Portion	5.0 mm
Working Length	380 mm
Minimum Instrument Channel Width	≥ 7.2 Fr
Additional Features	
Custom Buttons	PC200-AS, PC200-AR
Irrigation Tube Included	All Models

**PUSEN**

Zhuhai Pusen Medical Technology Co., Ltd.

**mediform**

MEDIFORM, spol. s r.o., Olšová 1, CZ Brno  
637 00 TELEFON: +420 541 221 392  
E-MAIL: [objednavky@mediform.cz](mailto:objednavky@mediform.cz) /  
[www.mediform.cz](http://www.mediform.cz)