

PUSEN


mediform

Single Use Flexible Video Ureteroscope

With Surgeon Controlled Pressure System

PU3033H | PU3033AH

*Precision Controlled Suction
Innovative Urological Solutions*



Integrated Suction Port

- Facilitate intrarenal pressure and temperature maintaining a safe threshold
- Optimal irrigation flow during laser activation
- Low-risk of postoperative infections and other complications

Optical System	
Field of View	120°
Direction of View	0°
Depth of Field	3-50 mm
Insertion Portion	
Controllable Portion	270° Up 270° Down
Maximum Insertion Portion Width	2.7 mm
Diameter of Insertion Portion	7.5 Fr
Distal End Diameter	7.5 Fr
Working Length	650 mm
Minimum Instrument Channel Width	3.6 Fr
Additional Features	
Custom Buttons	PU3033H <input type="radio"/> PU3033AH <input checked="" type="radio"/>
Irrigation Tube Included	PU3033H <input checked="" type="radio"/> PU3033AH <input checked="" type="radio"/>




Straight Channel



Integrated T-Connector



Ergonomic Handle



Kink Resistant

mediform

MEDIFORM, spol. s r.o., Olšová 1, CZ Brno
637 00 TELEFON: +420 541 221 392
E-MAIL: objednavky@mediform.cz /
www.mediform.cz

PUSEN

Zuhai Pusen Medical Technology Co., Ltd.